Description

Cytori® Celase® GMP is a blended proteolytic enzyme designed for *in vitro* dissociation of adipose tissue from a variety of species. Celase GMP is a mixture of highly purified collagenase and neutral protease enzymes, formulated for efficient, gentle, and reproducible dissociation of nucleated cells from their surrounding tissue matrix. The target substrates for Celase GMP are the collagen and non-collagen proteins that comprise extracellular matrix and entrap nucleated cells. Celase GMP is aseptically processed, sterile filtered, and highly purified to assure the lowest levels of impurities. Celase GMP is provided within a clear glass vial which is contained within a foam-protected canister.

Formulation

Celase GMP is a white lyophilized consisting of aseptically filled, blended enzymes, and a small quantity of buffer salts. Celase GMP is not manufactured with mammalian tissues or other components.

Reconstitution Procedure

NOTE: Do not use bacteriostatic water for injection. This type of water contains preservatives that inhibit collagenase enzyme activity.

NOTE: Reconstitute the entire vial. Do not weigh individual aliquots of the lyophilized. The introduction of moisture into the vial results in a decline in enzymatic activity.

- Reconstitute the lyophilized Celase GMP blended enzyme into stock solution as follows:
 - Remove the plastic cap and swab the rubber septum with an alcohol wipe.
 - Using a sterile needle and syringe, pierce the rubber septum and add 5 mL of sterile Lactated Ringer's solution to the Celase GMP vial.
 - Invert the vial at least 10 times, making sure the enzyme is fully dissolved.
- Celase GMP should be used immediately after reconstitution and must be used within 4 hours or refrozen up to 2 times over 6 months when stored at -25 °C to -15 °C.

NOTE: Please refer to Cytori's species-specific tissue dissociation protocols for detailed instructions on the application and use of the reconstituted Celase GMP enzyme.

Sterility

Celase GMP is aseptically processed and sterile filtered.

Storage and Handling

Immediately store Celase GMP in a -25 °C to -15 °C freezer upon arrival, protected from extreme temperature conditions. The lyophilized enzyme is stable when stored unopened at -25 °C to -15 °C through the expiration date printed on the label. The reconstituted enzyme may undergo up to 2 freeze-thaw cycles over 6 months when stored at -25 °C to -15 °C.

Warnings

- Do not introduce air into the vial. This will create positive pressure which may result in fluid leakage through the septum.
- May be harmful if inhaled or swallowed.
- Possible skin/eye irritant in the lyophilized form.
- Reduce the Celase GMP enzyme by 100-Fold (minimum) after adipose tissue digestion to prevent deleterious effects to cells and/or tissues.
- For adipose tissue dissociation, use a temperature range of 35 °C to 39 °C. Although lower temperatures may be used, Celase GMP enzyme activity and the rate of adipose tissue dissociation will be slowed.
- Exposure of the Celase GMP enzyme to divalent cation chelators removes zinc and calcium, rendering the enzyme inactive.
- Do not use Celase GMP with modifying factors, such as serum, Bovine serum albumin or Fetal calf serum, or protease inhibitors.
- Do not re-sterilize.
- Degradation of reconstituted Celase GMP enzyme may occur following each freeze-thaw cycle performed after initial use. Degraded enzyme may result in under-performance of the adipose tissue digestion process and may lead to inadequate and/or improper cell preparation.
- Any attempts to re-sterilize Celase GMP enzyme may result in significant damage to the enzyme and/or loss of enzymatic activity. Damage to the enzyme would likely lead to inadequate adipose tissue digestion resulting in improper cell preparation.
- Discard and DO NOT USE previously opened or damaged Celase GMP enzyme, and use only enzyme that is packaged in unopened and undamaged containers.
- DO NOT USE if there is a potential loss of sterility of the Celase GMP enzyme.

Precautions

- Check the product appearance before use.
- Piercing the septum with needles having a diameter larger than 16 G may result in loss of septum integrity.
- Aseptic techniques must be used to minimize the possibility of contamination to the cell preparation at all times.
- Dissociation time, enzyme ratios, digestion temperature, and enzyme concentration all affect the tissue-dissociation outcome.
- Do not use Celase GMP for dissociation of tissue other than adipose tissue.
- Do not use Celase GMP in combination with other enzymes.
- Do not use Celase GMP in combination with: bacteriostatic water, chelating agents, or serum products. These items may inactivate the enzyme.

Possible Adverse Effects

Excessive Celase GMP enzyme contact with the tissue may have a deleterious effect on cell viability.

For Research Use Only

When Celase GMP is intended for cellular and tissue based product applications, the procedure must meet the applicable local laws and regulations, including the involvement of ethical committees, if necessary. Verification of product use for a specific application is the sole responsibility of the applicant.

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