



*Celbrush<sup>®</sup>*

*For Precision Delivery of Micro Droplets*



MANUFACTURER  
Cytori Therapeutics, Inc.  
3020 Callan Road  
San Diego, CA 92121, USA  
Tel: +1.858.458.0900  
info@cytori.com

CUSTOMER SERVICE  
USA: +1.858.875.5245  
EU: +44.1244.360.426  
AP: +81.3.6870.7500  
customerservice@cytori.com

RM-098-LIT-US\_B-0615

Cytori, the Cytori Logo, and Celbrush are trademarks or registered trademarks of Cytori Therapeutics, Inc. in the United States and other select countries.

© 2015 Cytori Therapeutics, Inc. All rights reserved.

## PRODUCT DESCRIPTION

Cytori's Celbrush is a stainless steel surgical instrument for precise delivery of micro droplets. The thumb-brush design gives the Celbrush a mechanical advantage that minimizes the build-up of pressure, and provides superior tactile feedback during tissue dispersion. This allows the surgeon, with a simple brush of his or her thumb, to exactly control the micro droplet tissue volume injected into the defect area without ever looking at the syringe, while simultaneously avoiding transient syringe obstruction and over-filling of the defect space.

The 10 mL Celbrush is designed to deliver approximately 0.50 mL of tissue for each full brush of the operator's thumb. For even more precision, the 3 mL Celbrush delivers approximately 0.15 mL per brush.

### KEY FEATURES AND BENEFITS

- Precision delivery of Micro droplets
- Quick loading
- Precise control of injected volume
- Minimizes clogging and over-filling
- 10 mL Celbrush minimizes surgeon hand fatigue
- 3 mL Celbrush for more precise injections



## FREQUENTLY ASKED QUESTIONS

### Do I need to buy a specific type of syringe for my Celbrush?

Yes, not all syringes fit the Celbrush. For our customer's convenience, Cytori provides a Celbrush Starter Kit, which includes one box (100/box) of compatible syringes, and an injection cannula for optimal use with the Celbrush. Using other brands of syringes is not recommended.

### Can I use any cannula with my Celbrush?

Yes, you can use any cannula with the Celbrush and not void the Celbrush warranty.

### Does the Celbrush come in different sizes?

Yes. There are 2 standard sizes, the 10 mL and the 3 mL Celbrush. The 10 mL Celbrush is designed to deliver approximately 0.50 mL of tissue for each full brush of the operator's thumb. For even more precision, the 3 mL Celbrush delivers approximately 0.15 mL per brush.

### How do I sterilize my Celbrush?

The Celbrush is a standard reusable and autoclavable surgical instrument. It is provided non-sterile and must be properly cleaned and sterilized prior to each use. Please see the instructions for use included with your Celbrush for more information regarding product care and autoclave instructions.

### How do I reorder syringes and cannulas?

You may purchase syringes or cannulas from Cytori by contacting Customer Service at +1.858.875.5245 or directly from the manufacturer.

## ORDERING INFORMATION

The 3 mL and the 10 mL Celbrush are each supplied in kits that include the assembled Celbrush device, a box of syringes specified to fit the Celbrush, and an injection cannula.

To order your Celbrush Starter Kit, email Customer Support at [customerservice@cytori.com](mailto:customerservice@cytori.com), call +1.858.875.5245, or visit us online at [cytori.com](http://cytori.com).

CATALOG NO.	DESCRIPTION
510/CRS	Celbrush 10 mL
503/CRS	Celbrush 3 mL
520/CRS	Celbrush 10 mL Starter Kit
521/CRS	Celbrush 3 mL Starter Kit
522/CRS	Celbrush 3 mL + 10 mL Starter Kit